

ALKYL PHENOLIC RESIN  
DPL.FP.SPER14

DRT- 4001

Sr.No	Test/Parameter	Specification	Method Reference
1	Appearance	Violet Coloured Lumps	DPL/CPR/001
2	Softening Point ( Ball and Ring)	90°C – 100°C	DPL/CPR/005
3	Methylol Content	12 – 16 %	DPL/CPR/011

**Application:** Used for making Neoprene/ Rubber based adhesives for shoe Industry.

ALKYL PHENOLIC RESIN  
DPL.FP.SPER15

DRT- 4002

Sr.No	Test/Parameter	Specification	Method Reference
1	Appearance	Pale Yellow Coloured Lumps	DPL/CPR/001
2	Softening Point ( Ball and Ring)	90°C – 100°C	DPL/CPR/005
3	Methylol Content	12 – 16 %	DPL/CPR/011

**Application:** Used for making Neoprene/ Rubber based adhesives for shoe Industry.